

ThermaStop™-RT

Improved RT-PCR Sensitivity, Specificity, and Product Yield




Figure 1. ThermaStop-RT greatly increases probe signals for detecting Hepatitis C Virus (HCV) RNA sequence using either two-step or one-step RT-PCR

ThermaStop-RT is a simple-to-use reagent for One-Step or Two-Step Reverse Transcription-PCR (RT-PCR) initiated with gene-specific primers. ThermaStop-RT interacts with reverse transcriptase at low temperatures to reduce priming errors that lead to non-specific products. By eliminating these errors, ThermaStop-RT significantly increases the yield of the desired specific products and increases assay sensitivity.

How to Use:

One unit of ThermaStop-RT is defined as the amount required for optimal performance in RT-PCR samples containing 50 units of reverse transcriptase and one unit of hot start Taq polymerase in a volume of 20 µl and a reverse transcription temperature of 50°C.

Reduced Non-Specific Amplification




Figure 2. ThermaStop-RT dramatically reduces non-specific amplification in one-step RT-PCR. Five samples without and five samples with ThermaStop-RT show (left to right) amplification of the mitochondrial nad5 mRNA from coconut palm leaves, synthetic coconut palm viroid RNA, or co-amplification of both nad5 mRNA and one of three serial dilutions of the viroid RNA. Only the specific 152 and 127 base pair products are detected in samples with ThermaStop-RT. Samples without the reagent have lower levels of those products and high levels of non-specific product. Bands on the left are from a quantitative DNA ladder (100-500 base pairs).

For Two-Step RT-PCR

Typical 20 µL cDNA synthesis mixture:

Reagent	Final Concentration	Volume/sample
RNase-free water		12.15 µl
5X RT buffer	1X	4 µl
10 mM dNTPs	0.4 mM	0.8 µl
5 U/µL ThermaStop-RT	0.05 U/µl	0.2 µl
200 U/µl RTase	2.5 U/µl	0.25 µl
10 µM primers	0.3 µM	0.6 µl
RNA template (e.g.)		2 µl

- Mix the first 5 components before adding primers and RNA template. Primers and RNA template can be prepared in a separate mixture with the RT buffer and incubated at an appropriate annealing temperature, if desired, before mixing with other components.
- Incubate cDNA synthesis mixture at 50°C or higher for 15-30 minutes, or 45°C for 30-60 minutes, then 95°C for 2 minutes.
- Dilute aliquots of the cDNA sample 4- to 10-fold with PCR reagent mixture and perform thermal cycling as usual.

For One-Step RT-PCR

Typical 20 μ L sample:

Reagent	Final Concentration	Volume/sample
RNase-free water		11.95 μ l (adjust as needed for 20 μ l final volume)
5X RT-PCR buffer	1X	4 μ l
10 mM dNTPs	0.4 mM	0.8 μ l
5 U/ μ L ThermaStop-RT	0.05 U/ μ l	0.2 μ l
200 U/ μ l RTase	2.5 U/ μ l	0.25 μ l
5 U/ μ l Hot-Start Taq	0.05 U/ μ l	0.2 μ l
10 μ M primers	0.3 μ M	0.6 μ l
RNA template		2 μ l (0.5 to 5 μ l range)

- It may be necessary to test a range of ThermaStop-RT concentrations (e.g. 0.5 to 4 units) for RT-PCR kits containing unknown Units of reverse transcriptase. ThermaGo should not be used with one-step RT-PCR.
- Mix the first 6 components before adding primers and RNA template. Primers and RNA template may be prepared in a separate mixture with the RT-PCR buffer and incubated at an appropriate annealing temperature, if desired, before mixing with other components.
- Incubate mixture at 50°C or higher for 15-30 minutes, or 45°C for 30-60 minutes, then 95°C for 2 minutes and a thermal cycling protocol appropriate for the primers and buffer conditions.

Applications of Interest:

ThermaStop-RT improves reverse transcription of virtually any RNA target. It is particularly useful in detecting low-copy targets or in assays that include multiple primer pairs (multiplexing). ThermaStop-RT has been tested with reverse transcriptases derived from MMLV and is likely to work with other RNA-dependent DNA polymerases. Please contact ThermaGenix Technical Support with questions regarding specific enzymes.

How to Prepare for Use:

ThermaStop-RT is shipped as a dry reagent in 250 or 1250 unit quantities and does not contain magnesium, dNTPs, or buffer components. Centrifuge tube briefly to insure the dried reagent is at the bottom of the tube. Add sterile, molecular grade 10 mM Tris, pH 8.3 in a volume to provide a final stock concentration of 5 Units/ μ l. Vortex 1-2 minutes, then centrifuge briefly. Allow tube to sit at room temperature for 15 minutes with occasional mixing to ensure reagent is completely dissolved.

Recommended Storage:

Stock solution should be refrigerated (4°C) or frozen (-20°C). Aliquot into smaller volumes as necessary to avoid freezing and thawing more than 5 times.

Results to Expect:

ThermaStop-RT will generate higher yields of the specific products (e.g. stronger probe signals, specific bands on gels) and less non-specific products compared to samples without ThermaStop-RT. The reagent should also improve detection of RNA targets at a much lower concentration.

Ordering Information

Part Number	Description	Number of standard reactions
rtTS-125	ThermaStop- RT 250 Units	250
rtTS-625	ThermaStop- RT 1250 Units	1250



16 Tech Circle Suite 201
Natick, MA 01760
(833) NEXT PCR
(507) NEXT PCR